MEDICONNECT

MEDICAL LOW PRESSURE HOSE LINE SYSTEMS

INSTRUCTIONS FOR USE

Document Number:	735 800 000 XXX
Date of Issue:	01/02/2008
Revision Number:	1

EN
DE
FR
NL
DA
FI
NO
SV
ES
PT
IT
CS
HR
HU
RU
EL
TR
ZH
JA
SK
PL

1. Foreword

GCE low pressure medical hoses are classified as Class IIa medical devices according to the Medical Device Directive 93/42/EEC. Their compliance with the requirements of the 93/42/EEC Medical Device Directive is based upon the EN 739 standard.

2. Intended use

Medical gas hoses are intended for the transfer of medical gases or vacuum, and/or the connection of two medical products (e.g. a pressure regulator with an emergency ventilator).

3. Operational, transport and storage safety requirements

	Operations	Transport	Storage
 <ul style="list-style-type: none"> • Keep the product and its associated equipment away from <ul style="list-style-type: none"> - heat sources (fire, cigarettes, ...), - flammable materials, - oil or grease, - water, - dust. 	✓	✓	✓
 <ul style="list-style-type: none"> • Always maintain oxygen cleanliness standards. 	✓	✓	✓
 <ul style="list-style-type: none"> • Use the product and its associated equipment in a well ventilated area. 	✓		

Statutory laws, rules and regulations for medical gases, accident prevention, and environmental protection must be observed.

4. Personnel instructions and training

According to Medical Devices Directive 93/42/EEC, the product owner must ensure that all personnel handling the product are provided with this "Instructions for Use" manual and are trained on product safety and use.

5. Connection of the hose line

First connect the hose to the device supplied. Connect the other end of the hose to the source.

Probe: ¹⁾(e.g. DIN 13260, NF, REN, etc.) Verify that the plug and coupling are not damaged and are marked with the same gas type code. When the plug is inserted, it engages audibly in the coupling.

NIST/DISS-connector: ¹⁾Verify that the NIST/DISS nipple and the NIST/DISS housing are marked with the same gas type code, and that the screw joint and o-ring are not damaged.

- 
The nut must be hand tightened only! Do not use tools.

For more information regarding gas specific medical quick connectors, see Appendix 2.

Regularly inspect the hose during use to ensure that the hose is not torn, bent, twisted, or affected by excessive pull.

¹⁾ Depending on type

Use a clean, moistened cloth for cleaning the product. If necessary, disinfect the product with a commercial substance (e.g. isopropyl alcohol).

- **Please ensure that no disinfectant enters the hose, as there is a risk of respiratory injuries.**

7. Maintenance

1. Periodic checks:

Using a leak-detection liquid designated for medical purposes, check the tightness of the hose connections a minimum of twice a year. If a leak is detected, the entire hose line must be replaced. Also, twice a year verify that the hose has not exceeded the maximum use as explained in section 8 of this manual.

2. Inspection before use:

Before each use, perform an audible and visual check for leaks. If a leak is detected, or if there is damage to the connections or hose, cracks at the hose line, abrasion at the connecting parts (i.e. damaged or worn chrome), the entire hose line must be replaced.

Follow these maintenance procedures and local regulatory instructions, (e.g. France circulation letter 963921, dated 22nd July 1996).

3. Spare parts

In case of damage to the o-ring, the o-ring may be replaced with the following spare parts, as applicable:

O-ring NIST-plug O2: 14115260

O-ring NIST-plug AIR/O2: 14115261

O-ring DISS-nipple O2: 14115262

8. Marking

Colour coding of the hose and, if applicable, the connectors is in accordance with EN739 or per customer's request. The hose is marked for the gas for which it is intended.

If the GCE marking is not visible on both sleeves of the hose, the hose is not original equipment and must be replaced.

The lot number is written on one of the hose sleeves and product label, as follows: **XXXXXXXX/MMYY**
XXXXXXXX - product number, MM - month of production, YY - year of production (2 last numbers of the year).

If the hose is in good condition, without any leakage or damage, the maximum use of the hose line system is 10 years (use last two numbers of lot number plus 10).

9. Other

Disposal of old or damaged products must comply with national regulatory.

10. Glossary

Consult operating instruction

Keep away from heat and flammable material

Keep dry

Suitable for Hospital care use

Caution

Keep away from oil and grease

Upper and lower humidity limit

Upper and lower temperature limit

Suitable for Emergency care use

Product part number

Date of manufacturing


Use by date

Appendix 2, Quick coupling feature and connecting/disconnecting procedure

Manufactured by
GCE s.r.o.
Zizkova 381
583 81 Chotebor
Czech Republic

Tel : +420 569 661 111
Fax : +420 569 661 602
<http://www.gce-online.com>
© GCE s.r.o.

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